

FREEZING KIT

W / Penicillin-streptomycin
W / Bovine serum albumin
W / HEPES buffer



www.futurebiomedia.com

A. USAGE:

- For cryopreservation of sperm for further use in Assisted Reproductive Technologies.

B. CONTENTS:

- Freezing solution (FS).
- Product insert.

C. STORAGE:

- The kit Keeps dry at **2 – 8 °C**.

D. STERILITY:

- Sterile Filtered.

E. COMPOSITION:

Inorganic salts

Sodium Chloride
Sodium Phosphate
Potassium Chloride
Magnesium Sulfate
Calcium Chloride

Buffers

Sodium Bicarbonate
HEPES

Amino Acids

Glycine

Antibiotics

Penicillin-streptomycin

Energy Sources

Dextrose

Protein

Bovine serum albumin

Cryoprotectant

Sucrose

Glycerol

Water

WFI Quality

F. PRECAUTIONS:

- Respect storage conditions of the product.
- Do not use the product after its expiry date.
- Manipulate the product in aseptic conditions (e.g.: under laminar air flow).
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall, etc.).

G. PROCEDURES:

- BioMEDIA provide easy cryopreservation protocol for sperms originated from ejaculate or TESE extracted in simplified way to obtain high survival rate **60 – 80 %**.

Sperm cryopreservation protocol:

- Cryopreservation is performed on native semen samples or processed samples.
- Ensure Freezing solution (FS) is well mixed at room temperature before use.
- For native semen sample allow the semen to liquefy at **37 °C** for **30 minutes**.
- Mix equal volume of semen with Freezing solution (FS).
- Add the Freezing solution (FS) in drops gently to avoid osmolarity shock within **2 minutes** with continues mixing and leave the mixture for **10 minutes** at room temperature for equilibration.

6. Suck the sample/Freezing solution (FS) mixture into the freezing straw (e.g. **SPERM CRYO SYSTEM**), leaving approximately **1.5 cm** of air at the end of the straw and seal the straw.
7. Place the straw horizontally in a liquid nitrogen vapor for freezing occur.
8. leave for **30 minutes**, transfer straws quickly into liquid nitrogen and store at **-196 °C**.

Sperm thawing protocol:

1. Remove straw from the liquid nitrogen and place the straw in tap water for **1 minutes (room temperature or 37 °C)**.
2. Cut the end of straw and place the open side inside a container and cut the head of straw for obtaining sample/Freezing solution (FS) mixture.
3. Dilute the sample/Freezing solution (FS) mixture in a suitable HEPES buffered media (e.g. **S-VIVO MEDIA**) at least 3 ml per **0.25 ml or 0.5 ml** of sample/Freezing solution (FS) mixture for removing cryoprotectant.
4. Centrifuge at **300 g** for **10 minutes**.
5. Resuspend pellet in a suitable volume of HEPES buffered media (e.g. **S-VIVO MEDIA**) according to pellet size.
6. Finally, asses recovered sample.

H. TECHNICAL ASPECTS:

1. Increasing of exposure time to Freezing solution (FS) may be harmful for sperm.
2. During thawing, removed the straw quickly from liquid nitrogen and put into tap water (**room temperature or 37 °C**) to avoid heat shock.
3. Quick dilution after thawing is very important to avoid toxicity and increase survival rate.

