




SCIENTIFIC ANALYSIS INSTRUMENTS


# Laser ToF LT2 PLUS



## Microbiology Milestone

# Revolutionizing Microbial Diagnostics

The LT2 Plus is a compact, benchtop MALDI-TOF mass spectrometry system engineered to streamline microbial identification in clinical, research, and industrial environments. In combination with the BactoSCREEN-ID software, it delivers fast, cost-effective, and highly accurate identification of Gram-positive and Gram-negative bacteria, yeasts, and molds.



## Rapid Identification

Identifies 96 samples in under 15 minutes, which is significantly faster than traditional methods.

Requires minimal user experience with a simple workflow: spot bacteria, add matrix, and analyze.

## Ease of Use



## High Accuracy

Provides identification accuracy exceeding 98%, ensuring reliable results for various applications.

Offers low consumable costs compared to traditional identification methods, saving resources.


## Cost-Effective



# Simplified Workflow

## User-Friendly, Efficient, and Scalable


### MALDI Analysis Workflow



Additional steps are available for yeast & mold identification.  
No prior experience is needed to operate the system.

# User Friendly Software

The BactoSCREEN-ID software is specifically created to quickly determine the taxonomic classification of a micro-organism being studied at the species level. This is achieved by comparing its mass spectrometry profile to the profiles of known strains of micro-organisms present in the corresponding database.



The results are presented in a clear manner that allows users to identify the species level along with its corresponding score value. Additionally, the color codes are assigned based on the score value to further enhance clarity.

## Score Value Meaning

Organism (Genus , Species)	Score Result
Staphylococcus Saprophyticus LT7AB4	0.91
Staphylococcus xylosus LT3593	0.57
Staphylococcus kloosii 7A72	0.48



High Confidence Species

Range from 0.79 to 1.0 indicates high confidence identification for species. This is represented by the color green.

Range from 0.49 to 0.79 indicates high confidence genus identification. Species identification has low confidence and is represented by yellow.


High Confidence Genus




No Identification Possible

Range from 0 to 0.49 indicates that no identification is possible. This is represented by the color red.

Dendrogram of the selected wells



With Bactoscreen, users can create their own database, conduct statistical analyses, generate dendrograms, and compare data.



# Rapid Blood Culture Workflow

## Make a Critical Difference for Your Sepsis Patients

The Rapid Workflow is designed to meet the urgent need for fast and reliable identification from Positive Blood Cultures (PBC).

This solution empowers clinical microbiologists to deliver results that can significantly impact patient care, particularly in time-sensitive cases like sepsis.

Following identification, the system allows for spectral analysis to detect potential resistance markers—offering early alerts about antimicrobial resistance without requiring additional procedures.



### Benefits of a Blood culture workflow



#### Rapid Identification

Quick species-level identification after PBC alert.



#### Faster Decisions

Physicians gain early, actionable results promptly.



#### Subculturing Bypass

Eliminates subculturing, saving diagnostic process time.




#### Efficient Preparation

PBC samples prepared quickly and efficiently.



#### Resistance Detection

Early detection of resistance for timely intervention.



# Accurate Identification of Molds and Multicellular Fungi

## Overcoming a Longstanding Challenge

Identifying filamentous fungi and molds has traditionally been a complex task in microbiology.

Even advanced MALDI-TOF systems face limitations, as fungal growth conditions can significantly impact spectral results, leading to inconsistent mass spectra. Through dedicated research and development, we have successfully generated reproducible reference spectra, enabling reliable and consistent identification of these challenging organisms.

### Preparing Fungi for Identification



# LT2 Plus Specifications

## Power in a Compact Design



### Mass Resolution

Resolution greater than 2,500 m/ $\Delta$ m (FWHM).



### Mass Range

Mass range extends up to 500,000 Da.



### Accuracy

Accuracy is less than 100 ppm (internal).



### Sensitivity

Sensitivity reaches below 1 fmol detection limit.



### Laser

Pulsed UV laser operating at 337 nm.



### Vacuum

Air-cooled TMP with oil-free diaphragm pump.



### Stage

Motorized stage designed for 96-well plates.



### Software

Software provides full MS control and database support.



### Dimensions

Physical dimensions are 52×83×130 cm.



### Power

Power consumption is less than 1 kVA.

# ACCESSORIES & CONSUMABLES



## HCCA MALDI MATRIX KIT

High Purity Matrix Substance enables easy and convenient preparation of MALDI matrix solution for MALDI-TOF- MS measurement of peptides and proteins. Each Kit contains : 4 vials for 1000 samples.



## MALDI AUTOMATIC SPOTTER

is designed for an automated fast and accurate sample preparation for microbial identification .The spotter can process up to two MALDI targets simultaneously, requiring only one second per sample spot. For higher throughput, the instrument has an automatic target drying function that provides a fast and efficient sample preparation.



## BACTERIA STANDARD KIT

The Bacteria Standard consists of E. coli DH5 proteins that have been spiked with two additional proteins in the upper mass range, allowing for calibration over a mass range of 3 to 18 kDa.

Each Kit contains : 4 vials for 200 Calibration .



SCIENTIFIC ANALYSIS INSTRUMENTS



## CONTACT US



+44 (0) 161 850 1231



[www.saiman.co.uk](http://www.saiman.co.uk)



Media House, Wright  
Old Trafford, Manchester  
M16 9HB, United Kingdom

