English

Quilted Heel Protector



9514

The **NYOrtho SUPER Quilted Heel Protector** is designed to cushion the heel, reduce pressure and protect against skin breakdown. It is made of quilted fabric with an extra-thick filling that redistributes pressure and minimizes friction and shear.

Adjustable hook-and-loop straps allow a secure, comfortable fit. The fabric is breathable, water resistant and antimicrobial, wicks moisture and allows airflow.

Models 9514 and 9514L, sold as a pair.

Soft hypoallergenic filling | **Padded closure** | **Water resistant** | **Machine washable** | **Antimicrobial**

SUPER Quilted Heel Protector | Models 9514, 9514L

Details

Product family: NYOrtho SUPER Quilted

Product name: Quilted heel protector (SUPER Quilted Heel Protector)

Models: 9514 (Standard) and 9514L (Bariatric), sold as a pair

Intended use: Recommended for patients with cognitive and/or physical disabilities that prevent them from performing protective behaviors such as repositioning or expressing discomfort

Features

- **SUPER FLUFFY filling:** soft, hypoallergenic microfiber for added comfort.
- **Quilted fabric with extra-thick filling** redistributes pressure and minimizes friction and shear.
- **Breathable, water-resistant, antimicrobial fabric** wicks moisture and allows airflow.
- **Padded hook-and-loop closure** allows easy application and removal.
- **Machine washable.**
- **Meets CA 117 standards.**

Specifications

Model 9514	Standard, 25 × 25 cm, sold as a pair
Model 9514L	Bariatric, 30 × 30 cm, sold as a pair
Filling	Hypoallergenic microfiber
Fabric	Breathable, water resistant, antimicrobial
Closure	Padded hook & loop
Laundering	Machine wash, up to 160°F / 71°C
Standard	CA 117
Approval	Approved by the Israeli Ministry of Health

Application and care

- Before applying, inspect the heels and lower legs for pressure injuries, redness, skin breakdown or swelling. The patient should preferably be lying down or seated with legs elevated.
- The heel protector is not a wound dressing and must never be applied directly over an open or untreated wound. If a wound is present, follow facility wound care protocol to clean, treat and dress it before applying the heel protector.
- Find the care label to identify the top of the heel protector. Open the strap, place the heel protector under and around the foot so that the heel rests on the internal cushioning surface, and fasten the strap across the top of the foot securely enough to keep the device in place but not so tight that circulation is impaired.
- Reassess the fit and the skin condition under the device at least every 2 hours or per facility protocol, and remove it periodically.
- Machine wash warm (up to 160°F / 71°C), dry on low heat, bleach as directed. Inspect after each wash and replace if the cushioning becomes compressed or the fabric, seams or straps are compromised.

Manufacturer: NYOrtho, USA
Importer: Medtex (Gur Prym Agencies Ltd.)
11 Gush Etzion St., Givat Shmuel, Israel
Tel: 050-7213213 · info@medtex.co.il · www.medtex.co.il

REF08112026

rev. 261108