



PRODUCT
BROCHURE

State-of-the-art aesthetic
equipment backed by expertise.
Trusted by professionals, tested
for reliability, and designed for
success.

LUMA
SOLUTIONS




About Luma Solutions

At Luma Solutions, we specialize in providing high-quality aesthetic machines designed to ensure safety and client satisfaction. With experience of over three decades in business and extensive expertise in the aesthetic industry, our team's journey began with the challenge of finding reliable equipment for our own clinics. This challenge inspired us to conduct thorough research on the best available options.

Our Commitment to Quality

As clinic owners ourselves, we understand the importance of utilizing effective machinery. We not only design a selection of our machines in the UK, but we also use several of them in our own practices. We are committed to selling only those machines that have a high client demand for treatments, ensuring that you can maximize your profits while providing sought-after services.



Trusted Manufacturers & Engineering Expertise

To guarantee excellence in our products, we meticulously choose our manufacturers based on their extensive quality control measures.

Their comprehensive teams of quality check engineers and dedicated safety departments ensure that every machine meets the highest standards before it ever reaches your clinic. Additionally, our engineers conduct thorough inspections and safety checks of all machines at our UK depot to verify their performance and safety prior to dispatching them to you. This rigorous attention to detail gives us confidence that the equipment we provide is both reliable and safe.

Through extensive discussions with certification bodies in the UK and collaboration with skilled, UK-based laser engineers who have worked with leading UK laser companies, we identified crucial criteria for selecting reputable manufacturers. These criteria include valid certifications, strong reputations, and a solid safety check department.

Our engineers now partner with Luma Solutions to offer maintenance checks, servicing, and repairs, ensuring your machines operate at peak performance.

Custom-Tailored Solutions

Our team works closely with these manufacturers, supplying detailed specifications and guidelines that allow us to tailor our machines precisely to the needs of our clients, offering unique solutions that stand out in the market

Reliable and Affordable Machines

We are proud to present a carefully sourced selection of top quality machines that meet industry standards. Our mission is to provide you with affordable and reliable machines that enhance your practice and drive profitability.

Comprehensive Training & Support

Additionally, Luma Solutions collaborates with experienced educators who possess a wealth of knowledge in the operation of aesthetic devices. They offer comprehensive, fully accredited device training, either at your premises or online. This ensures that you and your staff are fully equipped to utilize the technology effectively

Customizable payment solutions

Select from Various Budget friendly options to enhance your business's revenue. At Luma Solutions we provide tailored payment plans that align with your cash flow needs Enabling you to offer exceptional skin care treatments with cutting edge beauty devices .

Available payment options

- One -time purchase: Make a 50% deposit to secure your order then pay the remaining balance upon delivery.
- Simple rental: easily rent a premium device and connect with us at the end of your rental for the next steps .
- Flexible financing: explore various payment plans with our trusted finance partners.



LUMA

HYDRALUX FACIAL

Our non-invasive treatment improves skin health, addresses individual concerns, and creates a glow like no other. Unlike traditional facials, it's performed via device rather than by hand, for deeper penetration and longer lasting results.



TECHNICAL SPECIFICATIONS:

Voltage	100-120V/ 60Hz 200-240V/50Hz
Display	10.4 inches
Power	500W
RF Frequency	0.2-0.5MHz
Heat&Cool	Heat: 36-46 Cool: 0-6
Ultrasound	1Mhz
Aqua Peel Probes	Soft: 3 pcs Hard: 8 pcs
Aqua Peel	Vacuum intensity: 4-5kg Vacuum Pressure: 96Kpa Water flow: ≤50L
Solution Bottle	A: Cleaning B: Whitening and moisturizing C: Deep layer cleaning D: Pipe cleaning
CLEANSING SHOVEL Frequency	26-30MHz
EMS NANO Mesotherapy Size	4X4mmX250Um
MICRO LIFT	Frequency: 200Hz Current: 50-800Ua
Skin Analyzer	28 MILLION HD MEGA PIXELS



Get the Best Skin
of Your Life



Meet all skin care needs!

- Fine lines and wrinkles
- Elasticity and firmness
- Even tone and vibrancy
- Skin texture
- Brown spots
- Oily and congested skin
- Enlarged pores

3 STEPS. 30 MINUTES.

Get the Best Skin of Your Life

It delivers instant and long-lasting results that you can see and feel.

STEP1 CLEANSE+PEEL

Uncover a new layer of skin with gentle exfoliation and relaxing resurfacing.

STEP2 EXTRACT+HYDRATE

Remove debris from pores with painless suction. Nourish with intense moisturizers that quench skin.


STEP3 FUSE+PROTECT

Saturate the skin's surface with antioxidants and peptides to maximize your glow.

DEEP CLEANSING SKIN SYSTEM

Cleaning blackheads, acne, and mites while moisturizing nutrients supplying; fining pores, increasing the skin transparency, and improving skin texture.

What is special about it?
Three times the power of the operating head.




The three major areas with clear division of work can achieve three times the effect at one time.

ANTI-ALLERGY

SOOTHING

REPAIR

SKIN TIGHTENING



INTELLIGENT OPERATING SYSTEM

- ⦿ Versatile Functions
- ⦿ Versatile Languages
- ⦿ User - Friendly UI
- ⦿ 10.4" display
- ⦿ 10 in 1 Platform

VERSATILE APPLICATORS - 10IN 1PLATFORM

1.AQUA PEEL



Use vacuum to open pores for a deep cleanse and easy extraction of impurities and dead skin cells. It also allows the active ingredients in our boosters to penetrate deeper, for the most effective treatment possible.

2.OXYGEN JET



High-pressure oxygen injection, small particles and strong penetration Penetrates deep into the basal layer, enhances water retention, and absorbs more easily by narrowing pores

3.EMS FACE



Nanocrystals promote the absorption of skin care products by opening channels for skin absorption Therefore, when nanocrystals are used, skin care products can be more absorbed by the skin

4.COLD-HEAT



Quickly lower the temperature below 0 degrees for skin quality, and achieve the effect of removing redness, anti-sensitivity, firming and shaping

5.MICRO LIFT



Anion and cation EMS microcurrent alternates between positive and negative, so that cations are imported and anions are exported improve facial contours

6.OXYGEN



Create an oxygen supply environment, activate cells, and enhance absorption Inject the product into the deep layer of the skin and cooperate with fruit acid skin peeling products

7.ULTRASOUND



Ultrasound—1 to 3 million vibrations per minute, opens the pores, allowing the essence to quickly penetrate into the base layer of the skin, allowing the absorption rate to reach more than 90%

8.RF



The probe can emit radio waves that vibrate 6 million times per second. The high-frequency vibrations are the molecules in the skin that generate heat due to friction, so as to achieve the effect of firming the skin and improving wrinkles

9.CLEANSING SHOVEL



Use high-frequency vibration waves to soften the stratum corneum and deeply clean the skin dirt

10.ANALYZER



Skin Analyzer:
The best tool to make skin analysis and offer best solutions !

A Type


- 1.Aqua Peel
- 2.A Type(OXYGEN JET)
- 3.EMS FACE
- 4.COLD-HEAT
- 5.MICRO LIFT
- 6.OXYGEN
- 7.ULTRASOUND
- 8.RF
- 9.CLEANSING SHOVEL
- 10.ANALYZER

B Type

- 1.Aqua Peel
- 2.B Type(OXYGEN JET)
- 4.COLD-HEAT
- 5.MICRO LIFT
- 6.OXYGEN
- 7.ULTRASOUND
- 8.RF
- 9.CLEANSING SHOVEL


HOW EMS FACE WORKS?

Nanocrystals use nano-scale contacts, which are as thin as more than 250 nanometers and as short as 0.1 mm, which can just open the stratum corneum of the skin without touching the dermis. It not only solves the problem of skin absorption, but also achieves painless and non-invasive treatment!




Natural application and absorption

In the case of natural application and absorption of skin care products, the molecules of skin care products only penetrate into the epidermis. It is difficult to penetrate into the deeper skin layer, and the absorption rate is very low.



Nanocrystalline nanochip (The product can penetrate evenly into the dermis.)

Nanochips only penetrates the barrier layer of the epidermis without touching blood vessels and nerves. It is painless and non-invasive, and the absorption capacity is greatly improved.



EMS function

The microelectronic heat conduction technology opens the pores, which can increase the product absorption rate by more than 40%.

Principle of nanochip introduction


- ⦿ Nanocrystals are made of monocrystalline silicon. In addition to being tough, monocrystalline silicon can also be thin enough and short enough to only break through the barrier layer of epidermal absorption, and it's far from touching blood vessels and nerves, so you can't feel pain, nor will it bleed.
- ⦿ Nanocrystals promote the absorption of skin care products by opening the skin absorption channel, so skin care products can be more absorbed by the skin when nano crystals are introduced.




Sagging skin



Dry and water shortage



Fine water wrinkles



Gloomy





LIFE-CHANGING RESULTS

The amazing results are both instant and long-lasting.

FINE LINES



before after

NASOLABIAL FOLDS




before after

BROWN SPOTS



before after

OILY & CONGESTED SKIN



before after

ENLARGED PORES



before after

EVENTONE



before after

SAGGING SKIN



before after

DRY SKIN



before after