

# POWER HAPE 2<sup>™</sup>

Multi-functional Body Shaping System

Bi-polar & Multi-polar RF, Vacuum, and low level laser

#### **SPECIFICATION**

Radio Frequency	2MHz(Bi-polar, Multi-polar, Thermal)
Laser Spectrum	658nm
Vacuum Range	280mmHg
Applicators	Roll Handpiece, Multi-Polar Handpiece, Thermal Handpiece
Input Power	AC 100~240V, 50/60Hz
Panel	10.2 inch Touch LCD
Dimension(WxDxHmm)	414 X 523 X 1,230



Head office: 120, Gieopdosi-ro, Jijeong-myeon, Wonju-si, Gangwon-do, Korea.

Tel: +82-33-760-1700(Rep.) Fax: +82-33-746-9716

Seoul office: #1701~1702, 145, Gasan digital 1-ro, Geumcheon-gu, Seoul, Korea

Tel: +82-2-514-9713(Rep.) Fax: +82-2-514-9716

E-mail: es@esglobal.co.kr / www.esglobal.co.kr





PowerShape 2 is an integrated system for body contouring and face lifting combined with Multi & Bi-polar RF, Thermal monopolar RF, Vacuum and low level laser



## Optimum Body Manager

Four functions in one system

### > Radio Frequency

Radio Frequency delivers thermal energy deep into dermal layer, and involves applying electromagnetic waves to return physiological processes.

It offers effective treatment results in cellulite, body sculpting and collagen regeneration by various RF technologies of bipolar and multipolar RF


#### Vacuum


Increases blood circulation and evens out the skin to optimize efficient energy delivery.

Massaging effect facilitates lymphatic drainage and reshapes the body smoother.

#### > Low Level Laser

Delivers specific Wavelength of light(658nm) which liquefies the fat for easier fat removal and significantly reduces the size of fat cells.







\*Reference: Plastic And reconstructive Surgery, September 1, 2002. Fat liquefaction: Effect of low level laser energy on Adipose Tissue by Neira et al.

Non-Surgical No Pain No Side Effect No-Discomfort

### > Applicators



Large size roll handpiece pulls the target area skin to enhance vacuum suction effect





Multi-Polar RF technology painlessly delivers focused and optimal RF energy to the skin surface and to subcutaneous fat layers.









By delivering Thermal energy in Non-ablative method, RF energy penetrates deep into the dermis layer, resulting in collagen stimulation without any damage to the superficial skin.



### > Strengths

Absolute treatment result

Various applicators

- Powerful movement by Dual Motorized Rollers

4 types of vibration pulse

10.2 inch wide screen