

**F R E E D O M<sup>®</sup>**




**K N E E**



**Freedom Total Knee<sup>®</sup> System  
Partial Knee**



**Pursue Life<sup>™</sup>**



### KEEL CUTS

- Improve implant strength and stability providing a mechanical lock
- Provides even cement coverage throughout the underside of implant

### CEMENT DAM

- Cement containment to minimize extrusion

## THE LATEST INNOVATION IN PARTIAL KNEE TECHNOLOGY

The Freedom<sup>®</sup> Knee Partial Knee provides the following:

- Faster Recovery Time
- Minimized Post Operative Pain
- Early Mobilization
- Minimal Bone resection
- Easy, reproducible procedure with guided resection technique. Little to no free hand technique
- Full femoral resurfacing



Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee<sup>®</sup> System.

**CAUTION:** THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

**WARNING:** THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Freedom Knee<sup>®</sup> is manufactured by Maxx Orthopedics, Inc. FREEDOM and FREEDOM KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

© 2016 Maxx Orthopedics. All rights reserved. Updated February 2016.



Maxx Orthopedics, Inc.  
531 Plymouth Road, Ste. 526  
Plymouth Meeting, PA 19462 USA

Confidential. For illustration only. Material not yet approved for distribution.